
CT-Scan protocol - requirements for the CT data set

- Scan format: DICOM file format
- Scan direction: Axial - *performed in maximal inspiration*
- Scan window/filtering: Tissue
- Reconstructions: no reconstruction required – preferably the raw data set
- Slice thickness: ≤ 1.2 mm
- Slice distance: 0.6 mm (equal distance)

File labeling

(For reasons of data protection law, we can only accept data records in anonymized form)
Please rename the DICOM dataset according following standard

YYYYMMDD_Prescriber_123456789_01
Date_Name of Prescriber_Pat.ID or Acronym_Consecutive no

Data transfer options – please specify

- We want to upload dataset (After receiving the request Novatech SA need to send an upload link)

- We send CD-ROM via air mail to: **Novatech SA**, Customer Service, Gustav-Krone-Strasse 7, 14167 Berlin, Germany

- We pass CD-ROM to Mr./Ms.:

Company:

Please send this form to Email: intl-sales@bessgroup.com **Fax:** +49 30 816 909 16
